

BioCytex

140, Chemin de l'Armée d'Afrique
13010 Marseille
France

Tel : +33 (0) 4 96 12 20 40
Fax : +33 (0) 4 91 47 24 71
Email : info@biocytex.fr

Sept.06

DATA SHEET

CD4-FITC (Cat # 5104-F100T)

MOUSE MONOCLONAL ANTIBODY ANTI-HUMAN CD4, FITC CONJUGATED

Clone	EDU.2
Isotype	Mouse IgG2a, Kappa.
Partner of fusion	Mouse myeloma NS1 cells.
Immunogen	Stimulated human leukocytes.
Specificity	The antibody recognizes the CD4 antigen of 55 kD, presents on T-helper/inducer lymphocytes and monocytes.
Application	Flow cytometry.
Form	FITC conjugated purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 2 mL.
Size	100 tests, ready for use.
Suggested amount	20 µL / test per 10 ⁶ cells in 100 µL sample.
Preservative	Sodium azide < 0,1%.
Storage	The conjugated antibody should be stored in the dark at +2-8°C. Do not freeze.
HLDA Workshop	Leukocyte Typing II, vol (1), Springer-Verlag (eds), pp 16, 25 (1986).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC PROCEDURE