

BioCytex

140, Chemin de l'Armée d'Afrique
13010 Marseille
France
Tel : +33 (0) 4 96 12 20 40
Fax : +33 (0) 4 91 47 24 71
Email : info@biocytex.fr

Mar 05

DATA SHEET

CD62P-FITC

(Cat # 5111-F100T)

MOUSE MONOCLONAL ANTIBODY ANTI HUMAN P-SELECTIN (CD62P), FITC CONJUGATED

Clone	LYP20
Isotype	IgG1, Kappa
Partner of fusion	Sp2/0-Ag14
Immunogen	Chymotrypsin-treated platelets
Specificity	This antibody recognizes human P-selectin also known as CD62P and GMP140. It has been also described to cross-react with rat P-selectin ⁽²⁻³⁾ . CD62P is expressed on platelets, megakaryocytes and endothelial cells.
Application	Flow cytometry
Form	FITC conjugated purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 2 mL.
Size	100 tests, ready for use.
Suggested amount	20 µL/test per 10 ⁶ cells in 100 µL sample.
Preservative	Sodium azide < 0.1%.
Storage	The conjugated antibody should be stored in the dark at +2-8°C. Do not freeze.
References	<ol style="list-style-type: none">1. Puente Navazo MD. <i>et al.</i>, Thromb Haemost. 2001, vol 86, 2, 680-685.2. Ruchaud-Sparagano MH. <i>et al.</i>, Biochem J. 1998, 332, 309-314.3. Chignier E. <i>et al.</i>, Thromb Haemost. 1994, vol 72, 5, 745-749.4. Parmentier S. <i>et al.</i>, Blood 1991, vol 77, 8, 1734-1739.
Licence	Non exclusive licence with INSERM U331

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURE