

BioCytex

140, Chemin de l'Armée d'Afrique
13010 Marseille
France
Tel : +33 (0) 4 96 12 20 40
Fax : +33 (0) 4 91 47 24 71
Email : info@biocytex.fr

Aug. 03

DATA SHEET

Mouse IgG1-PE

(Cat #5108-PE100T)

MOUSE IgG1 ISOTYPIC CONTROL, R-PE CONJUGATE

Clone	2DNP 2H11/2H12
Isotype	IgG1 Kappa
Partner of fusion	SP2/0-Ag14
Immunogen	DNP-BSA
Specificity	DNP (2,4-Dinitrophenyl) is an antigen not expressed on normal human leukocytes, platelets and red cells.
Application	Flow cytometry.
Form	R-Phycoerythrin-conjugated purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 2 mL.
Size	100 tests, ready for use.
Suggested amount	20 µL/test for 100 µL of sample.
Preservative	Sodium azide < 0.1%.
Storage	The conjugated antibody should be stored in the dark at +2-8°C. Do not freeze.
HLDA Workshop	N/A

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURE