

BioCytex

140, Chemin de l'Armée d'Afrique
13010 Marseille
France
Tel : +33 (0) 4 96 12 20 40
Fax : +33 (0) 4 91 47 24 71
Email : info@biocytex.fr

Nov. 05

DATA SHEET

Mouse IgG2b (Cat #5117-P)

MOUSE IgG2b ISOTYPIC CONTROL, PURIFIED

| | |
|--------------------------|--|
| Clone | 2DNP 14G5 |
| Isotype | IgG2b, kappa |
| Partner of fusion | P3X63Ag8/653 |
| Immunogen | DNP-BSA |
| Specificity | This antibody recognizes DNP (2,4-Dinitrophenyl). This molecule is not expressed on normal human leukocytes, platelets and red cells. |
| Application | Flow cytometry. |
| Form | Purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 1 mL. |
| Size | 0.1 mg |
| Suggested amount | Use this isotypic control at the same concentration than the specific monoclonal antibody tested. |
| Preservative | Sodium Azide < 0.1%. |
| Storage | As delivered at +2-8°C until expiration date. After opening, aliquote contents and freeze at -20°C. Avoid repeated freezing and thawing. |
| HLDA Workshop | N/A |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURE.