

BioCytex

140, Chemin de l'Armée d'Afrique
13010 Marseille
France
Tel : +33 (0) 4 96 12 20 40
Fax : +33 (0) 4 91 47 24 71
Email : info@biocytex.fr

Mar 05

DATA SHEET

CD42b-FITC

(Cat# 5113-F100T)

MOUSE MONOCLONAL ANTIBODY ANTI HUMAN CD42b, FITC CONJUGATED

Clone	ALMA 19
Isotype	IgG1, Kappa
Partner of fusion	P3X63 Ag8.653
Immunogen	Human platelets
Specificity	This antibody recognizes CD42b, platelet glycoprotein Ib alpha (GpIb α). CD42b is expressed on platelets and megakaryocytes.
Application	Flow cytometry
Form	FITC conjugated purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 2 mL.
Size	100 tests, ready for use.
Suggested amount	20 μ L/test per 10^6 cells in 100 μ L sample.
Preservative	Sodium azide < 0.1%.
Storage	The conjugated antibody should be stored in the dark at +2-8°C. Do not freeze.
Reference	Perrault C. <i>et al.</i> , Blood 2003, 101, 3477-3484.
Licence	Non exclusive licence with the "Etablissement Français du Sang" - Alsace

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURE