

# BioCytex

140, Chemin de l'Armée d'Afrique  
13010 Marseille  
France  
Tel : +33 (0) 4 96 12 20 40  
Fax : +33 (0) 4 91 47 24 71  
Email : info@biocytex.fr

Mar 05

## DATA SHEET

# Mouse IgG1

(Cat # 5108-P)

### MOUSE IgG1 ISOTYPIC CONTROL, PURIFIED

<b>Clone</b>	2DNP 2H11/2H12
<b>Isotype</b>	IgG1, Kappa
<b>Partner of fusion</b>	SP2/0-Ag14
<b>Immunogen</b>	DNP-BSA
<b>Specificity</b>	DNP (2,4-Dinitrophenyl) is an antigen not expressed on normal human leukocytes, platelets and red cells.
<b>Application</b>	Flow cytometry.
<b>Form</b>	Purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 1 mL.
<b>Size</b>	0.1 mg/mL
<b>Suggested amount</b>	Use this isotypic control at the same concentration than the specific monoclonal antibody tested.
<b>Preservative</b>	Sodium azide < 0.1%.
<b>Storage</b>	As delivered at +2-8°C until expiration date. After opening, aliquote contents and freeze at -20°C. Avoid repeated freezing and thawing.
<b>HLDA Workshop</b>	N/A

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURE