

BioCytex

140, Chemin de l'Armée d'Afrique
13010 Marseille
France
Tel : +33 (0) 4 96 12 20 40
Fax : +33 (0) 4 91 47 24 71
Email : info@biocytex.fr

Nov. 05

DATA SHEET

Mouse IgG2a-PE

(Cat# 5127-PE100T)

MOUSE IgG2a ISOTYPIC CONTROL, R-PE CONJUGATED

Clone	2DNP 16C12
Isotype	IgG2a, Lambda
Partner of fusion	P3X63Ag8/653
Immunogen	DNP-BSA
Specificity	This antibody recognizes DNP (2,4-Dinitrophenyl). This molecule is not expressed on normal human leukocytes, platelets and red cells.
Application	Flow cytometry.
Form	R-Phycoerythrin conjugated purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 2 mL.
Size	100 tests, ready to use.
Suggested amount	20 μ L/test per 10^6 cells in 100 μ L sample.
Preservative	Sodium azide < 0.1%.
Storage	The conjugated antibody should be stored in the dark at +2-8°C. Do not freeze.
HLDA Workshop	N/A

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURE.