

BioCytex

140, Chemin de l'Armée d'Afrique
13010 Marseille
France
Tel : +33 (0) 4 96 12 20 40
Fax : +33 (0) 4 91 47 24 71
Email : info@biocytex.fr

May 05

DATA SHEET

Mouse IgG2a

(Cat #5127-P)

MOUSE IgG2a ISOTYPIC CONTROL, PURIFIED

Clone	2DNP16C12
Isotype	IgG2a, Lambda
Partner of fusion	SP2/0-Ag14
Immunogen	DNP-BSA
Specificity	DNP (2,4-Dinitrophenyl) is an antigen not expressed on normal human leukocytes, platelets and red cells.
Application	Flow Cytometry
Form	Purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 1 mL.
Size	0.1 mg
Suggested amount	Use this isotypic control at the same concentration than the specific monoclonal antibody tested. Quantitative flow cytometry : follow the package insert of CELLQUANT Calibrator kit (BioCytex ref. 7208).
Preservative	Sodium Azide < 0.1%
Storage	As delivered at +2-8°C until expiration date. After opening, aliquote contents and freeze at -20°C. Avoid repeated freezing and thawing.
HLDA Workshop	N/A

FOR RESEARCH USE ONLY NOT FOR USE IN DIAGNOSTIC PROCEDURE