

BioCytex

140, Chemin de l'Armée d'Afrique
13010 Marseille
France
Tel : +33 (0) 4 96 12 20 40
Fax : +33 (0) 4 91 47 24 71
Email : info@biocytex.fr

Nov. 05

DATA SHEET

Mouse IgG2b-FITC

(Cat #5117-F100T)

MOUSE IgG2b ISOTYPIC CONTROL, FITC CONJUGATED

Clone	2DNP 14G5
Isotype	IgG2b, kappa
Partner of fusion	P3X63Ag8/653
Immunogen	DNP-BSA
Specificity	This antibody recognizes DNP (2,4-Dinitrophenyl). This molecule is not expressed on normal human leukocytes, platelets and red cells.
Application	Flow cytometry.
Form	FITC conjugated purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 2 mL.
Size	100 tests, ready to use.
Suggested amount	20 µL/test per 10 ⁶ cells in 100 µL sample.
Preservative	Sodium azide < 0.1%.
Storage	The conjugated antibody should be stored in the dark at +2-8°C. Do not freeze.
HLDA Workshop	N/A

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURE.