

# BioCytex

140, Chemin de l'Armée d'Afrique  
13010 Marseille  
France  
Tel : +33 (0) 4 96 12 20 40  
Fax : +33 (0) 4 91 47 24 71  
Email : [info@biocytex.fr](mailto:info@biocytex.fr)

Nov. 05

## DATA SHEET

# Mouse IgG2b-PE

(Cat #5117-PE100T)

MOUSE IgG2b ISOTYPIC CONTROL, R-PE CONJUGATED

<b>Clone</b>	2DNP 14G5
<b>Isotype</b>	IgG2b, kappa
<b>Partner of fusion</b>	P3X63Ag8/653
<b>Immunogen</b>	DNP-BSA
<b>Specificity</b>	This antibody recognizes DNP (2,4-Dinitrophenyl). This molecule is not expressed on normal human leukocytes, platelets and red cells.
<b>Application</b>	Flow cytometry.
<b>Form</b>	R-Phycoerythrin conjugated purified immunoglobulin in PBS-BSA 0.1%, pH 7.2, liquid, 2 mL.
<b>Size</b>	100 tests, ready to use.
<b>Suggested amount</b>	20 µL/test per 10 <sup>6</sup> cells in 100 µL sample.
<b>Preservative</b>	Sodium azide < 0.1%.
<b>Storage</b>	The conjugated antibody should be stored in the dark at +2-8°C. Do not freeze.
<b>HLDA Workshop</b>	N/A

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURE.